

CIPAC

COLLABORATIVE INTERNATIONAL PESTICIDES ANALYTICAL COUNCIL LIMITED


Commission Internationale des Méthodes d'Analyse des Pesticides (CIMAP)

CHAIRMAN:

Dr. Ralf Hänel
Federal Office of Consumer Protection and Food Safety
Messeweg 11/12
38104 Braunschweig, Germany
Telephone: +49 531 299 3506
Telefax: +49 531 299 3002
E-Mail: ralf.haenel@bvl.bund.de

TREASURER:

Brian E. Hocken
17 Claygate Avenue
Harpenden, Herts, AL5 2HE, England
Telephone: +44 1582 764277
Telefax: +44 1582 461625
E-Mail: brian@hockene.freemove.co.uk



SECRETARY:

László Bura
European Food Safety Authority
Largo N. Palli 5/A, I-43121 Parma, Italy
Telephone: +39 0521 036827
Telefax: +39 0521 0360827
E-Mail: laszlo.bura@efsa.europa.eu

ASSISTANT SECRETARY:

Sonia Tessier
Chemicals Regulation Directorate of HSE
Mallard House, Kings Pool, 3,
Peasholme Green
YO 17PX York, UK
Tel +44 1904455905
E-Mail: sonia.tessier@hse.gsi.gov.uk

Internet: <http://www.cipac.org>

CIPAC INFORMATION SHEET No 285

Your Ref:
Our Ref: LB
15 September 2010

578. flumioxazin

This is to inform you that JAPAC is planning to organize a collaborative trial on a method for flumioxazin.

Principle. The content of flumioxazin is determined by reversed phase high performance liquid chromatography on a 250 x 4.6 mm ODS column (Phenomenex Gemini 5µ C18 (5 µm), or equivalent) with mobile phase composed of acetonitrile/water, 50/50 (v/v) and UV detection at 288 nm. Quantification is done by external standardization.

Scope. It is intended to analyze five samples consisting of three technical materials and two WG formulations. Samples and flumioxazin analytical standard will be supplied.

Time schedule. The collaborative trial is intended to be carried out during November 2010 – February 2011.

Those who wish to participate should contact the organizer by 30th October 2010.

Contact:

Mr. Yasushi Asada
Sumitomo Chemical Co., Ltd.
Organic Synthesis Research Laboratory
Analytical Science Group
3-1-98 Kasugadenaka Konohanaku,
Osaka, 554-8558 Japan

Tel: +81-6-6466-5171
Fax: +81-6-6466-5427
E-mail: asada@sc.sumitomo-chem.co.jp

Yours sincerely

L. Bura
Secretary